


FINESSE INJECTABLE™



SUMMA

The Revascularization Multi-Tool

The Finesse Injectable™ PTA Balloon Dilatation Catheter is the only available device that combines a state-of-the-art low-profile 0.014" platform angioplasty balloon with on-the-fly injection capability while over-the-wire, allowing Finesse to function as a crossing catheter, diagnostic catheter, and angioplasty balloon catheter--all without catheter exchanges.

Technical Specifications

Catheter type	OTW	Balloon lengths	20-225 mm
Recommended guidewire	0.014"	Balloon type	Semi-compliant
Tip	0.016" entry profile	Distal Shaft	2.9F (0.038")
Balloon material	Pebax	Injection ports	3 ports proximal to balloon
Balloon coating	Hydrophilic	Usable length	65, 100, 150 cm
Balloon diameters	2.0--4.0 mm	Introducer sheath	5F

Ordering Information

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200201503	00810017490004	2	20	150	5	10	12	0.014
FS250201503	00810017490011	2.5	20	150	5	8	12	0.014
FS300201503	00810017490028	3	20	150	5	9	12	0.014
FS350201503	00810017490035	3.5	20	150	5	8	12	0.014
FS400201503	00810017490042	4	20	150	5	7	12	0.014
FS200401503	00810017490059	2	40	150	5	10	12	0.014
FS250401503	00810017490066	2.5	40	150	5	8	12	0.014
FS300401503	00810017490073	3	40	150	5	9	12	0.014
FS350401503	00810017490080	3.5	40	150	5	8	12	0.014
FS400401503	00810017490097	4	40	150	5	7	12	0.014
FS200601503	00810017490103	2	60	150	5	10	12	0.014
FS250601503	00810017490110	2.5	60	150	5	8	12	0.014
FS300601503	00810017490127	3	60	150	5	9	12	0.014
FS350601503	00810017490134	3.5	60	150	5	8	12	0.014
FS400601503	00810017490141	4	60	150	5	7	12	0.014
FS201001503	00810017490158	2	100	150	5	10	12	0.014
FS251001503	00810017490165	2.5	100	150	5	8	12	0.014
FS301001503	00810017490172	3	100	150	5	9	12	0.014
FS351001503	00810017490189	3.5	100	150	5	8	12	0.014
FS401001503	00810017490196	4	100	150	5	7	12	0.014
FS202251503	00810017490202	2	225	150	5	10	12	0.014
FS252251503	00810017490219	2.5	225	150	5	8	12	0.014
FS302251503	00810017490226	3	225	150	5	9	12	0.014
FS352251503	00810017490233	3.5	225	150	5	8	12	0.014
FS402251503	00810017490240	4	225	150	5	7	12	0.014

For Healthcare Professionals only. Important information: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions.

Copyright Summa Therapeutics 2023. All rights reserved. The Summa Therapeutics logo and Finesse logo are trademarks of Summa Therapeutics and may be registered in the U.S. and other countries. All other marks are the property of their respective owners.

Global Headquarters
 Summa Therapeutics, LLC
 20 University Road, Ste. 500
 Cambridge, MA 02138
 Tel. (877) 267-2685
 info@summatherapeutics.com

www.summatherapeutics.com



Ordering Information

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200201003	00810017490257	2	20	100	5	10	12	0.014
FS250201003	00810017490264	2.5	20	100	5	8	12	0.014
FS300201003	00810017490271	3	20	100	5	9	12	0.014
FS350201003	00810017490288	3.5	20	100	5	8	12	0.014
FS400201003	00810017490295	4	20	100	5	7	12	0.014
FS200401003	00810017490301	2	40	100	5	10	12	0.014
FS250401003	00810017490318	2.5	40	100	5	8	12	0.014
FS300401003	00810017490325	3	40	100	5	9	12	0.014
FS350401003	00810017490332	3.5	40	100	5	8	12	0.014
FS400401003	00810017490349	4	40	100	5	7	12	0.014
FS200601003	00810017490356	2	60	100	5	10	12	0.014
FS250601003	00810017490363	2.5	60	100	5	8	12	0.014
FS300601003	00810017490370	3	60	100	5	9	12	0.014
FS350601003	00810017490387	3.5	60	100	5	8	12	0.014
FS400601003	00810017490394	4	60	100	5	7	12	0.014
FS201001003	00810017490400	2	100	100	5	10	12	0.014
FS251001003	00810017490417	2.5	100	100	5	8	12	0.014
FS301001003	00810017490424	3	100	100	5	9	12	0.014
FS351001003	00810017490431	3.5	100	100	5	8	12	0.014
FS401001003	00810017490448	4	100	100	5	7	12	0.014
FS202251003	00810017490455	2	225	100	5	10	12	0.014
FS252251003	00810017490462	2.5	225	100	5	8	12	0.014
FS302251003	00810017490479	3	225	100	5	9	12	0.014
FS352251003	00810017490486	3.5	225	100	5	8	12	0.014
FS402251003	00810017490493	4	225	100	5	7	12	0.014

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200200653	00810017490509	2	20	65	5	10	12	0.014
FS250200653	00810017490516	2.5	20	65	5	8	12	0.014
FS300200653	00810017490523	3	20	65	5	9	12	0.014
FS350200653	00810017490530	3.5	20	65	5	8	12	0.014
FS400200653	00810017490547	4	20	65	5	7	12	0.014
FS200400653	00810017490554	2	40	65	5	10	12	0.014
FS250400653	00810017490561	2.5	40	65	5	8	12	0.014
FS300400653	00810017490578	3	40	65	5	9	12	0.014
FS350400653	00810017490585	3.5	40	65	5	8	12	0.014
FS400400653	00810017490592	4	40	65	5	7	12	0.014
FS200600653	00810017490608	2	60	65	5	10	12	0.014
FS250600653	00810017490615	2.5	60	65	5	8	12	0.014
FS300600653	00810017490622	3	60	65	5	9	12	0.014
FS350600653	00810017490639	3.5	60	65	5	8	12	0.014
FS400600653	00810017490646	4	60	65	5	7	12	0.014
FS201000653	00810017490653	2	100	65	5	10	12	0.014
FS251000653	00810017490660	2.5	100	65	5	8	12	0.014
FS301000653	00810017490677	3	100	65	5	9	12	0.014
FS351000653	00810017490684	3.5	100	65	5	8	12	0.014
FS401000653	00810017490691	4	100	65	5	7	12	0.014
FS202250653	00810017490707	2	225	65	5	10	12	0.014
FS252250653	00810017490714	2.5	225	65	5	8	12	0.014
FS302250653	00810017490721	3	225	65	5	9	12	0.014
FS352250653	00810017490738	3.5	225	65	5	8	12	0.014
FS402250653	00810017490745	4	225	65	5	7	12	0.014

Global Headquarters
Summa Therapeutics, LLC
20 University Road, Ste. 500
Cambridge, MA 02138
Tel. (877) 267-2685
info@summaterapeutics.com

www.summaterapeutics.com

