


# FINESSE BTK MULTICATHETER™



**SUMMA**

# The Revascularization Multi-Tool

The Finesse BTK Multicatheter is the only available device that functions as a state-of-the-art low-profile crossing catheter, an injectable diagnostic catheter, and a high-performance small vessel angioplasty balloon catheter, all-in-one, all while maintaining guidewire access.

## Technical Specifications

Catheter type	OTW	Balloon lengths	20-225 mm
Recommended guidewire	0.014"	Balloon type	Semi-compliant
Tip	0.016" entry profile	Distal Shaft	2.9F (0.038")
Balloon material	Pebax	Injection ports	3 ports proximal to balloon
Balloon coating	Hydrophilic	Usable length	65, 100, 150 cm
Balloon diameters	2.0--4.0 mm	Introducer sheath	5F

## Ordering Information

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200201502	00810017490752	2	20	150	5	10	12	0.014
FS250201502	00810017490769	2.5	20	150	5	8	12	0.014
FS300201502	00810017490776	3	20	150	5	9	12	0.014
FS350201502	00810017490783	3.5	20	150	5	8	12	0.014
FS400201502	00810017490790	4	20	150	5	7	12	0.014
FS200401502	00810017490806	2	40	150	5	10	12	0.014
FS250401502	00810017490813	2.5	40	150	5	8	12	0.014
FS300401502	00810017490820	3	40	150	5	9	12	0.014
FS350401502	00810017490837	3.5	40	150	5	8	12	0.014
FS400401502	00810017490844	4	40	150	5	7	12	0.014
FS200601502	00810017490851	2	60	150	5	10	12	0.014
FS250601502	00810017490868	2.5	60	150	5	8	12	0.014
FS300601502	00810017490875	3	60	150	5	9	12	0.014
FS350601502	00810017490882	3.5	60	150	5	8	12	0.014
FS400601502	00810017490899	4	60	150	5	7	12	0.014
FS201001502	00810017490905	2	100	150	5	10	12	0.014
FS251001502	00810017490912	2.5	100	150	5	8	12	0.014
FS301001502	00810017490929	3	100	150	5	9	12	0.014
FS351001502	00810017490936	3.5	100	150	5	8	12	0.014
FS401001502	00810017490943	4	100	150	5	7	12	0.014
FS202251502	00810017490950	2	225	150	5	10	12	0.014
FS252251502	00810017490967	2.5	225	150	5	8	12	0.014
FS302251502	00810017490974	3	225	150	5	9	12	0.014
FS352251502	00810017490981	3.5	225	150	5	8	12	0.014
FS402251502	00810017490998	4	225	150	5	7	12	0.014

For Healthcare Professionals only. Important information: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions.

Copyright Summa Therapeutics 2023. All rights reserved. The Summa Therapeutics logo and Finesse logo are trademarks of Summa Therapeutics and may be registered in the U.S. and other countries. All other marks are the property of their respective owners.

Global Headquarters  
Summa Therapeutics, LLC  
20 University Road, Ste. 500  
Cambridge, MA 02138  
Tel. (877) 267-2685  
info@summatherapeutics.com

[www.summatherapeutics.com](http://www.summatherapeutics.com)

**SUMMA**  
Rethink Angioplasty™

# Ordering Information

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200201002	00810017490257	2	20	100	5	10	12	0.014
FS250201002	00810017490264	2.5	20	100	5	8	12	0.014
FS300201002	00810017490271	3	20	100	5	9	12	0.014
FS350201002	00810017490288	3.5	20	100	5	8	12	0.014
FS400201002	00810017490295	4	20	100	5	7	12	0.014
FS200401002	00810017490301	2	40	100	5	10	12	0.014
FS250401002	00810017490318	2.5	40	100	5	8	12	0.014
FS300401002	00810017490325	3	40	100	5	9	12	0.014
FS350401002	00810017490332	3.5	40	100	5	8	12	0.014
FS400401002	00810017490349	4	40	100	5	7	12	0.014
FS200601002	00810017490356	2	60	100	5	10	12	0.014
FS250601002	00810017490363	2.5	60	100	5	8	12	0.014
FS300601002	00810017490370	3	60	100	5	9	12	0.014
FS350601002	00810017490387	3.5	60	100	5	8	12	0.014
FS400601002	00810017490394	4	60	100	5	7	12	0.014
FS201001002	00810017490400	2	100	100	5	10	12	0.014
FS251001002	00810017490417	2.5	100	100	5	8	12	0.014
FS301001002	00810017490424	3	100	100	5	9	12	0.014
FS351001002	00810017490431	3.5	100	100	5	8	12	0.014
FS401001002	00810017490448	4	100	100	5	7	12	0.014
FS202251002	00810017490455	2	225	100	5	10	12	0.014
FS252251002	00810017490462	2.5	225	100	5	8	12	0.014
FS302251002	00810017490479	3	225	100	5	9	12	0.014
FS352251002	00810017490486	3.5	225	100	5	8	12	0.014
FS402251002	00810017490493	4	225	100	5	7	12	0.014

Catalog Number	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Sheath Size (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Guidewire (in)
FS200200652	00810017490509	2	20	65	5	10	12	0.014
FS250200652	00810017490516	2.5	20	65	5	8	12	0.014
FS300200652	00810017490523	3	20	65	5	9	12	0.014
FS350200652	00810017490530	3.5	20	65	5	8	12	0.014
FS400200652	00810017490547	4	20	65	5	7	12	0.014
FS200400652	00810017490554	2	40	65	5	10	12	0.014
FS250400652	00810017490561	2.5	40	65	5	8	12	0.014
FS300400652	00810017490578	3	40	65	5	9	12	0.014
FS350400652	00810017490585	3.5	40	65	5	8	12	0.014
FS400400652	00810017490592	4	40	65	5	7	12	0.014
FS200600652	00810017490608	2	60	65	5	10	12	0.014
FS250600652	00810017490615	2.5	60	65	5	8	12	0.014
FS300600652	00810017490622	3	60	65	5	9	12	0.014
FS350600652	00810017490639	3.5	60	65	5	8	12	0.014
FS400600652	00810017490646	4	60	65	5	7	12	0.014
FS201000652	00810017490653	2	100	65	5	10	12	0.014
FS251000652	00810017490660	2.5	100	65	5	8	12	0.014
FS301000652	00810017490677	3	100	65	5	9	12	0.014
FS351000652	00810017490684	3.5	100	65	5	8	12	0.014
FS401000652	00810017490691	4	100	65	5	7	12	0.014
FS202250652	00810017490707	2	225	65	5	10	12	0.014
FS252250652	00810017490714	2.5	225	65	5	8	12	0.014
FS302250652	00810017490721	3	225	65	5	9	12	0.014
FS352250652	00810017490738	3.5	225	65	5	8	12	0.014
FS402250652	00810017490745	4	225	65	5	7	12	0.014

Global Headquarters  
Summa Therapeutics, LLC  
20 University Road, Ste. 500  
Cambridge, MA 02138  
Tel. (877) 267-2685  
info@summaterapeutics.com

[www.summaterapeutics.com](http://www.summaterapeutics.com)

